

RF QIMT

RF QAS

Innovation and Accuracy in CardioVascular Disease prevention

QIMT-QAS

CardioVascular Disease
and Prevention ▼

RF-data
technology ▼


RFQIMT ▼

RFQAS ▼


Bibliography

Introduction


- CardioVascular Disease is the number one cause of death and disability (*World Health Report 2004, WHO*)
- According to World Health Organization (WHO) estimates, in 2003, 16.7 million people around the globe die of CVD each year. This is over 29 percent of all deaths globally. (*www.who.int*)
- Each year CVD causes over 4.35 million deaths in Europe and over 1.9 million deaths in the EU. CVD causes nearly half of all deaths in Europe (49 percent) and in the EU (42 percent). (*European Cardiovascular Disease Statistics. 2005 Edition, British Heart Foundation Health Promotion Research Group*)
- In US, nearly 2300 Americans die of CVD each day, an average of 1 death every 38 seconds (*AHA Statistical Update Heart Disease and Stroke Statistics - 2010 Update. A report from the American Heart Association*)



CVD deaths
around the **globe**



CVD deaths
in **Europe**



CVD deaths
in **U.S.**

CardioVascular Disease and Prevention

Introduction

Why is prevention of CVD
needed?

Objectives
of CVD prevention

Guidelines
of CVD prevention

RF-data technology

RFQIMT

RFQAS

Bibliography

1/4

Why is prevention of CVD needed?

Because CVD is an important cause of:

- Premature mortality and morbidity
- Years of life lost (YLL) and Disability adjusted lifeyears (DALY's)
- Increasing health care costs



CardioVascular Disease and Prevention

Introduction

Why is prevention of CVD needed?

Objectives of CVD prevention

Guidelines of CVD prevention

RF-data technology

RFQIMT

RFQAS

Bibliography

2/4

Objectives of CVD prevention

Early detection reduces the risk of major Cardiovascular events resulting in reduced premature morbidity and mortality, prolongs good health and improves quality of life. QIMT and QAS can help assess the need for lifestyle changes, management of cardiovascular risk factors and the prophylactic use of certain drugs.

CardioVascular Disease and Prevention

[Introduction](#)[Why is prevention of CVD needed?](#)[Objectives of CVD prevention](#)[Guidelines of CVD prevention](#)[RF-data technology](#)[RFQIMT](#)[RFQAS](#)[Bibliography](#)

3/4



Contact Info





Home

Guidelines on CVD Prevention

The 2012 European Guidelines on CVD Prevention in Clinical Practice is a new shorter and more practical format around the 5 key questions of prevention:

- What is it?
- Why is it needed?
- For whom?
- By whom?
- How can be provided?

The European Guidelines state that **vascular ultrasound screening is reasonable for risk assessment in asymptomatic individuals at moderate risk, i.e. for primary prevention.**



CardioVascular Disease and Prevention

Introduction

Why is prevention of CVD needed?

Objectives of CVD prevention

Guidelines of CVD prevention

RF-data technology

RFQIMT

RFQAS

4/4

Bibliography


 Click to enlarge

Guidelines on CVD Prevention

The 2012 European Guidelines on CVD Prevention in Clinical Practice is a new shorter and more practical format around the 5 key questions of prevention:

- What is it?
- Why is it needed?
- For whom?
- By whom?
- How can be provided?

The European Guidelines state that **vascular ultrasound screening is reasonable for risk assessment in asymptomatic individuals at moderate risk, i.e. for primary prevention.**



CardioVascular Disease and Prevention

Introduction

Why is prevention of CVD needed?

Objectives of CVD prevention

Guidelines of CVD prevention

RF-data technology

RFQIMT

4/4

RFQAS

Bibliography




Click to enlarge

RF-technology

The RadioFrequency (RF) is received unfiltered from the ultrasound transducer; it is raw data and contains 100% of the reflected signal.

The Esaote Quality Arterial Stiffness (^{RF}QAS) and Quality Intima media Thickness (^{RF}QIMT) are based on RF-data technology. This allows an accurate assessment of arterial stiffness and vessel wall thickness.



While traveling through the body, ultrasonic waves react uniquely to different tissue types; the RF preserves all information content.

CardioVascular Disease
and Prevention ▼

RF-data
technology

RF-technology

RF-technology versus
video-based methods

Highlights
and advantages

^{RF}QIMT ▼

^{RF}QAS ▼


Bibliography

1/3


RF-technology versus video-based methods

The RF based technology provides several advantages over the old video-based method:



- The information available with radiofrequency contains more than double the information contained in the video signal
- The real-time character of the functional information has also an inherent measurement quality feedback by means of a B-mode overlay
- Due to the real-time feedback the operator is able to optimize probe position to achieve the best possible scan plane.



QIMT: 778 μ m



RF Signal



Assessment of Blood Vessel Wall Properties by means of ultrasound

Quality Intima Media

CardioVascular Disease
and Prevention

RF-data
technology

RF-technology

RF-technology versus
video-based methods

Highlights
and advantages

RFQIMT

RFQAS

2/3

Bibliography

Highlights and advantages

- Accuracy
- Quality feedback control
- Ease of use
- Decreased operator variability
- Real-time

CardioVascular Disease
and Prevention ▼

RF-data
technology

RF-technology

RF-technology versus
video-based methods

Highlights
and advantages

RFQIMT ▼

RFQAS ▼

3/3

Bibliography

Contact Info




Quality Intima Media Thickness


Intima-media thickness (IMT) of the inner arterial wall has been shown to correlate with the degree of atherosclerosis measured by autopsy that, in turn, has been found to correlate with atherosclerotic burden in other arterial beds¹.

Consequently, carotid IMT is considered a surrogate marker of subclinical atherosclerosis. Increased carotid IMT is associated with CV risk factors, prevalent CVD and coronary artery atherosclerosis. Additionally, it has been shown to predict future CV events and death, independently of CV risk factors²⁻³.

Structural and functional ultrasound imaging such as arterial stiffness (QAS) and blood vessel wall thickness (QIMT) are emerging as the most representative of arteriosclerosis and atherosclerosis, and are widely validated for monitoring CVD patients and as independent markers of cardiovascular risk⁴⁻⁸.



Atherosclerosis is a condition in which plaque builds up on the inner lining of arteries.



 [Click to enlarge](#)

CardioVascular Disease
and Prevention ▼

RF-data
technology ▼

RFQIMT

Quality Intima
Media Thickness

Measurement
of the RFQIMT

RFQIMT graph

RFQIMT report

1/4

RFQAS ▼



Bibliography

Quality Intima Media Thickness

Intima-media thickness (IMT) of the inner arterial wall has been shown to correlate with the degree of atherosclerosis measured by autopsy that, in turn, has been found to correlate with atherosclerotic burden in other arterial beds¹.

Consequently, carotid IMT is considered a surrogate marker of subclinical atherosclerosis. Increased carotid IMT is associated with CV risk factors, prevalent CVD and coronary artery atherosclerosis. Additionally, it has been shown to predict future CV events and death, independently of CV risk factors²⁻³.

Structural and functional ultrasound imaging such as arterial stiffness (QAS) and blood vessel wall thickness (QIMT) are emerging as the most representative of arteriosclerosis and atherosclerosis, and are widely validated for monitoring CVD patients and as independent markers of cardiovascular risk⁴⁻⁸.



 Click to enlarge

CardioVascular Disease
and Prevention ▼

RF-data
technology ▼

RFQIMT

Quality Intima
Media Thickness

Measurement
of the RFQIMT

RFQIMT graph

RFQIMT report

RFQAS ▼


Bibliography

1/4

Measurement of the ^{RF}QIMT

The ^{RF}QAS and ^{RF}QIMT measurements are performed at the Common Carotid Artery.

Mannheim consensus suggests measuring IMT and stiffness on the common carotid artery, where the accuracy and reproducibility of the measurements are superior.



Click to enlarge

CardioVascular Disease
and Prevention ▼

RF-data
technology ▼

^{RF}QIMT

Quality Intima
Media Thickness

Measurement
of the ^{RF}QIMT

^{RF}QIMT graph

^{RF}QIMT report

2/4

^{RF}QAS ▼

Bibliography

Contact Info




Home

Measurement of the ^{RF}QIMT

The ^{RF}QAS and ^{RF}QIMT measurements are performed at the Common Carotid Artery.

Mannheim consensus suggests measuring IMT and stiffness on the common carotid artery, where the accuracy and reproducibility of the measurements are superior.



Click to enlarge

CardioVascular Disease
and Prevention ▼

RF-data
technology ▼

^{RF}QIMT

Quality Intima
Media Thickness

Measurement
of the ^{RF}QIMT

^{RF}QIMT graph

^{RF}QIMT report

^{RF}QAS ▼

Bibliography

2/4

Contact Info



Home

^{RF}QIMT graph

QIMT is based on the Howard or MESA tables.


The Howard table represents the QIMT values in micrometers taken from the ARIC study over age.

The MESA table represents the QIMT values taken from the MultiEthnic Study of Atherosclerosis (MESA).

The measured value is displayed in the QIMT graph and enables the doctor and patient to compare that value with the normal values obtained in one of the two studies cited.

Age in years

Any point below the line indicates a low risk



CardioVascular Disease
and Prevention ▼

RF-data
technology ▼

^{RF}QIMT

Quality Intima
Media Thickness

Measurement
of the ^{RF}QIMT

^{RF}QIMT graph

^{RF}QIMT report

3/4

^{RF}QAS ▼

Bibliography

CARDIOVASCULAR REPORT

Date : 21-Jan-2013

Gender : Female


Age : 25 Years
Systolic BP : 120 mmHg
Diastolic BP : 80 mmHg

QUALITY INTIMA MEDIA THICKNESS

RIGHT

QIMT	:	425 μ m	Diameter	:	6.36 mm
QIMT SD	:	9 μ m	Diameter SD	:	0.03 mm

CALCULATIONS



Normal QIMT	:	380 μ m	Very High Very Low
Maximum QIMT (right)	:	425 μ m	
Expected Age	:	32 Y	
Risk of CVD based on highest QIMT	:	High	

QUALITY ARTERIAL STIFFNESS

RIGHT



Distension	:	622 μ m	Brachial Psys	:	120 mmHg
Distension SD	:	24 μ m	Brachial Pdia	:	80 mmHg
Diameter	:	6.57 mm			
Diameter SD	:	0.07 mm			

Quality Arterial Stiffness

The stiffening of arteries is the most important cause of increasing both systolic and pulse pressure as well as decreasing diastolic pressure in people of 40 years and above.

This leads to cardiovascular complications, such as left ventricular hypertrophy, aneurysm formation and rupture⁹.

Hence a noninvasive determination of arterial stiffness during a routine ultrasound examination of the carotid arteries is a powerful clinical tool in CVD prevention.



 Click to enlarge

CardioVascular Disease
and Prevention ▼

RF-data
technology ▼

RFQIMT ▼

RFQAS

Quality Arterial
Stiffness

Measurement
of the RFQAS

RFQAS report

1/3



Bibliography

Quality Arterial Stiffness

The stiffening of arteries is the most important cause of increasing both systolic and pulse pressure as well as decreasing diastolic pressure in people of 40 years and above.

This leads to cardiovascular complications, such as left ventricular hypertrophy, aneurysm formation and rupture⁹.

Hence a noninvasive determination of arterial stiffness during a routine ultrasound examination of the carotid arteries is a powerful clinical tool in CVD prevention.



 Click to enlarge

CardioVascular Disease
and Prevention ▼

RF-data
technology ▼

RFQIMT ▼

RFQAS

Quality Arterial
Stiffness

Measurement
of the RFQAS



RFQAS report

1/3

Bibliography

Measurement of the ^{RF}QAS

QAS provides a list of standard parameters calculated combining the measured values (Distension, Distension Waveform and Diameter) with the Brachial Blood Pressure measured externally by using an automated system or standard cuff.



Healthy 23-year-old subject (A) – Healthy 65-year-old subject (B)
The 2nd peak due to wave reflection is more evident in the younger subject.



Click to enlarge

CardioVascular Disease
and Prevention ▼

RF-data
technology ▼

^{RF}QIMT ▼

^{RF}QAS

Quality Arterial
Stiffness

Measurement
of the ^{RF}QAS


^{RF}QAS report

2/3


Bibliography

Measurement of the ^{RF}QAS

QAS provides a list of standard parameters calculated combining the measured values (Distension, Distension Waveform and Diameter) with the Brachial Blood Pressure measured externally by using an automated system or standard cuff.



Healthy 23-year-old subject (A) – Healthy 65-year-old subject (B)
The 2nd peak due to wave reflection is more evident in the younger subject.



Click to enlarge

CardioVascular Disease
and Prevention ▼

RF-data
technology ▼

^{RF}QIMT ▼

^{RF}QAS

Quality Arterial
Stiffness

Measurement
of the ^{RF}QAS

^{RF}QAS report

2/3

Bibliography

RFQAS report

PWV: When the heart contracts it generates a pulse or energy wave that travels through the circulation.


The speed of travel of this pulse wave (Pulse Wave Velocity or PWV) is related to the arterial stiffness of the arteries.

The stiffer the artery is, the higher the PWV will be.

The final report shows:

- Distension values
- Stiffness values
- Local pressure values

RIGHT CCA QAS			STIFFNESS		
DISTENSION:	413	µm	DC:	0.02	1/kPa
SD:	12	µm	CC:	0.78	mm³/kPa
DIAMETER:	7.14	mm	α	3.80	
SD:	0.28	mm	β	7.82	
BrP sys:	130.0	mmHg	PWV	6.97	m/s
BrP dia:	80.0	mmHg			



LOCAL PRESSURE		
LOC Psys:	125.7	mmHg
LOC Pdia:	80.0	mmHg
P(T1):	123.4	mmHg
AP:	2.8	mmHg
Alx:	-2.25	%
ICP:	39	mm
ED:	244	mm

CardioVascular Disease
and Prevention

RF-data
technology

RFQIMT

RFQAS

Quality Arterial
Stiffness

Measurement
of the RFQAS

RFQAS report

3/3

Bibliography

Bibliography

1. Lian Engelen, Isabel Ferreira, Coen D. Stehouwer, Pierre Boutouyrie, and Stephane Laurent, on behalf of the Reference Values for Arterial Measurements Collaboration, *European Heart Journal*, (2012) doi:10.1093/eurheartj/ehs380, **Reference intervals for common carotid intima-media thickness measured with echo-tracking: relation with risk factors.**
2. Collaboration, *European Heart Journal* (2010) 31, 2338–2350 doi:10.1093/eurheartj/ehq165 Corresponding author, Pierre Boutouyrie, Department of Pharmacology and INSERM U970, Hôpital Européen Georges Pompidou, Assistance Publique - Hôpitaux de Paris, Université Paris-, **Determinants of pulse wave velocity in healthy people and in the presence of cardiovascular risk factors: 'establishing normal and reference values' The Reference Values for Arterial Stiffness'**
3. Touboul PJ, Hennerici MG, Meairs S, Adams H, Amarenco P, Bornstein N, Csiba L, Desvarieux M, Ebrahim S, Hernandez Hernandez R, Jaff M, Kownator S, Naqvi T, Prati P, Rundek T, Sitzler M, Schminke U, Tardif JC, Taylor A, Vicaut E, Woo KS., Mannheim carotid intima-media thickness and plaque consensus (2004-2006-2011). **An update on behalf of the advisory board of the 3rd, 4th and 5th watching the risk symposia, at the 13th, 15th and 20th European Stroke Conferences.** Mannheim, Germany, 2004, Brussels, Belgium, 2006, and Hamburg, Germany, 2011.
4. Pignoli P, **Ultrasound B-mode imaging for arterial wall thickness measurement.** *Atherosclerosis Reviews*, (1984) 12:177-184
5. Howard G, Sharrett AR, Heiss G, Evans GW, Chambless LE, Riley WA, et al., **Carotid artery intimal-medial thickness distribution in general populations as evaluated by B-mode ultrasound.** *Stroke*. 1993;24:1297 - 304
6. Joseph F. Polak, MD, MPH, Moyses Szklo, MD, DrPH; Richard A. Kronmal, PhD; Gregory L. Burke, MD, MSc; Steven Shea, MD, MS; Anna E. H. Zavodni, MD; Daniel H. O'Leary, MD, **The Value of Carotid Artery Plaque and Intima-Media Thickness for Incident Cardiovascular Disease: The Multi-Ethnic Study of Atherosclerosis**
7. Stephane Laurent¹, John Cockcroft, Luc Van Bortel, Pierre Boutouyrie, Cristina Giannattasio, Daniel Hayoz, Bruno Pannier, Charalambos Vlachopoulos, Ian Wilkinson and Harry Struijker-Boudier, on behalf of the European Network for Non-invasive Investigation of Large Arteries, **Expert consensus document on arterial stiffness: methodological issues and clinical applications.**
8. James H. Stein, MD, FASE, Claudia E. Korcarz, DVM, RDCS, FASE, R. Todd Hurst, MD, Eva Lonn MD, MSc, FASE, Christopher B. Kendall, BS, RDCS, Emile R. Mohler, MD, Samer S. Najjar, MD, Christopher M. Rembold, MD, and Wendy S. Post, MD, MS, ASE CONSENSUS STATEMENT **Use of Carotid Ultrasound to Identify Subclinical Vascular Disease and Evaluate Cardiovascular Disease Risk: A Consensus Statement from the American Society of Echocardiography Carotid Intima-Media Thickness Task Force Endorsed by the Society for Vascular Medicine**
9. J. Polak, S. Person en G. Wei, **Segment-specific association of carotid intima-media thickness with cardiovascular risk factors,** *Stroke*, vol. 41, pp. 9-15, 2010.
10. M. Lorenz, S. Kegler von en H. Steinmetz, **Carotid intima-media thickening indicates a higher vascular risk across a wide range. Prospective data from the Carotid Atherosclerosis Progression Study (CAPS),** *Stroke*, vol. 37, pp. 87-92, 2006.
11. M. F. O'Rourke, J. A. Staessen en C. Vlachopoulos, **Clinical applications of arterial stiffness; definitions and reference values,** *American Journal of Hypertension*, vol. 15, nr. 5, pp. 426-444, 2002.

CardioVascular Disease
and Prevention ▼

RF-data
technology ▼

RFQIMT ▼

RFQAS ▼

Bibliography

1/1



ISO13485:2003 ISO 9001:2008

Esaote S.p.A.

Via di Caciolle, 15 50127 Florence, Italy, Tel. +39 055 4229 1, Fax +39 055 4229 208, international.sales@esaote.com

Via A. Siffredi, 58 16153 Genoa, Italy, Tel. +39 010 6547 1, Fax +39 010 6547 275, info@esaote.com

FRANCE

Esaote Medical

ZA du Bel Air 10, rue de Témara,
78105 Saint-Germain-en-Laye

Tel. +33 1 8204 8900, Fax +33 1 3061 7210

info@esaote.fr

BRASIL

Braslian Direct Office

Rua Tomas Carvalhal, 711

04006-001 São Paulo SP

Tel. +55 11 2589 0533 Fax +55 11 2589 0527

leonardo.pili@esaote.com.br

GERMANY

Esaote Biomedica Deutschland GmbH

Max-Planck-Straße 27a

50858 Köln

Tel. +49 2234 688 5600, Fax +49 2234 967 9628

info@esaote.de

ARGENTINA

Esaote Latinoamérica S.A.

San Martín 551, Cuerpo 'C', Piso 8, (C1004AAK)
Buenos Aires

Tel. +54 11 4326 1832, Fax: +54 11 4328 1245

info@esaote.com.ar

SPAIN

Esaote España S.A.

C/. Pont Reixat, 5

08960 Sant Just Desvern, Barcelona

Tel. +34 93 473 2090, Fax +34 93 473 2042

info@esaote.es

INDIA

Esaote Asia Pacific Diagnostic Private Limited

DLF IT Park, A - 44 & 45, Tower- C, Ground Floor,

Sector- 62, Noida, Uttar Pradesh, India

Pin Code: 201 301

Tel. +91 120 4732444, Fax +91 120 4750148

info@esaote.in

THE NETHERLANDS AND BELGIUM

Esaote Benelux B.V.

Philipsweg 1

6227 AJ Maastricht

Tel. +31 43 3824650, Fax +31 43 3824651

benelux@esaote.nl

HONG KONG AND FAR EAST

Esaote China Ltd

18/F, 135 Bonham Strand Trade Centre,

135 Bonham Strand, Sheung Wan, Hong Kong

Tel. +852 2545 8386, Fax +852 2543 3068

esaote@esaotechina.com

UK

Esaote UK

14, Cambridge Science Park

Milton Road, Cambridge, CB4 0FQ

Tel. + 44 1223 424499, Fax + 44 709 288 0231

infoUK@esaote.com

CHINA

Esaote Shenzhen Medical Equipment

Room 2608, Tower B

Beijing Global Trade Center

36 North Third Ring Road East,

Dongcheng District, 100013, Beijing

Tel. +86 010 58257766, Fax +86 010 52257760

NORTH AMERICA

Esaote North America

8000 Castleway Drive,

Indianapolis, IN 46250

Tel. +1 317 813 6000, Fax +1 317 813 6600

inquire@esaoteusa.com

RUSSIAN FEDERATION AND CIS

Esaote S.p.A.

18 Leningradsky prospekt

Off. 5 and 6, Moscow 125040

Tel. +7 495 232 0205, Fax +7 495 232 1833

esaotemoscw@yandex.ru