

Sleek compact design
for portable exam environments




Tilt-and-swivel monitor for optimal
viewing during procedures



Customizable control panel to personalize
your user interface



Dual touch screens
streamline exam workflow



A Clear Vision

Ultrasound reimaged

Born of a vision to deliver meaningful design innovations that benefit the user, the Acclarix AX8 features a host of design breakthroughs that make day-to-day operation easy, fast and intuitive. The result is an elegant simplicity where form and function meet at the tips of your fingers.

- Unparalleled image quality with one-key optimization for all exams and all modes
- Innovative gesture control user interface sets a new standard in workflow
- Tilt-and-swivel, HD monitor for optimal viewing
- Completely sealed control panel to aid in maintaining infection control

Transducers

C5-2Q  2 - 5 MHz Curved Linear Array	C5-2MQ  2 - 5 MHz 3D/4D Mechanical Curved Array	MC8-4Q  4 - 8 MHz Micro Convex Array	E8-4Q  4 - 8 MHz Endocavity Curved Array
L17-7HQ  7 - 17 MHz High Frequency Linear Array	L12-5Q  5 - 12 MHz Linear Array	L17-7SQ  7 - 17 MHz High Frequency Compact Linear Array	MC9-3TQ  3 - 9 MHz Micro Convex Array
P7-3Q  3 - 7 MHz Phased Array	P5-1Q  1 - 5 MHz Phased Array	E10-3HQ  3 - 10 MHz Endocavity Curved Array (200°)	E10-3BQ  3 - 10 MHz Endocavity Curved Array (200°)

General Imaging Ultrasound in a Whole New Light.

Acclarix[™] AX8
Diagnostic Ultrasound System



GESTURE-CONTROL USER INTERFACE

APPLICATION-BASED WORKFLOW

USER-FOCUSED DESIGN

UNPARALLELED VALUE

DEFINITIVE IMAGE QUALITY





Global Headquarters:
Edan Instruments, Inc. | 15 Jinhui Road, Pingshan District, Shenzhen
518122 P.R. China | +86.755.26898326 | www.edan.com | info@edan.com
U.S. and Canada inquiries:
EDAN Diagnostics, Inc. | 9918 Via Pasar, San Diego, CA 92126
+1.858.750.3066 | www.edandiagnostics.com | edan-info@edandiagnostics.com




© Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice.
No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.

ENG-US-AX8-GI-V1.1-20200819







L12-5Q
Nerve




L17-7HQ
Breast



L12-5Q
Elastography



P5-1Q
Cardiac



E10-3HQ
Uterus



C5-2MQ
3D/4D

It's a Tour de Force

Redefining innovation through value and performance

It's the dawning of a new day in the world of compact ultrasound—and it turns out you really can have it all. The remarkable Acclarix™ AX8 Compact Ultrasound System delivers a powerhouse combination of features in a remarkably small package—making it the ideal choice for tight environments. The Acclarix AX8 has been designed from the ground up with a relentless focus on delivering unexpected levels of innovation and performance at a price point that is equally surprising.

- Definitive Imaging
- Distinctive Design
- Intelligent Workflow
- Intrinsic Quality

Stunning Clarity

Definitive image quality

Sleek design and user-friendly features are just part of the story. The Acclarix AX8 is also designed to deliver amazing clarity and stunning image quality quickly—eliminating the need to tweak imaging controls.

- High fidelity, high-channel count architecture results in superb detail resolution, particularly at depth
- Tissue Adaptive Imaging (TAI) continuously and automatically optimizes imaging allowing more focus on the patient
 - In B-mode, TAI fine tunes multiple parameters to provide the best possible image quality
 - In Doppler, TAI automatically adjusts for flow state providing improved continuity, border detection and fill-in



The Virtue of Value

How can Edan deliver so much innovation AND so much value? By capitalizing on the experience and expertise of a truly global enterprise. Acclarix technology is the result of Silicon Valley innovation and engineering excellence combined with Chinese design and manufacturing proficiency. The Acclarix AX8 delivers unmatched value and performance across a broad range of applications.

- Abdomen
- Difficult-to-Image
- OB/GYN
- Small Parts
- Breast
- Nerve
- Musculoskeletal
- Intraoperative
- Prostate
- Vascular
- Cardiac Screening